
Network Setup for TPTF Products

Note: You must have purchased a network license from BioEx for these instructions to work. If you have TPTF Office please do not use these instructions. The correct instructions can be found at <http://www.bioexsystems.com/documents.htm>.

Server Install

1. Insert the BioEx CD ROM into your drive on the server.
2. Click **Start: Run**
3. Type **d:\setup** (d = your CD drive) and then click the **OK** button.
4. During the setup, install your Maker product onto a drive which is accessible to all machines on the network that are going to run the TPTF Product
5. **Do Not** select the **Network Client** check box.
6. Make sure that users have full rights on the TPTF server directory.
7. Start your product and register it.

Workstation Install

1. Insert the BioEx CD ROM into your drive on the workstation.
2. Click **Start: Run**
3. Type **d:\setup** (d = your CD drive) and then click the **OK** button.
4. Install the TPTF product onto a local drive.
5. During the setup, select the **Network Client** check box.*
6. On the following screen, select the TPTF directory on the server.
7. After installing start the program and register it.

Enabling More than One Concurrent Copy on the Network

1. Start the TPTF product on any machine.
2. Open the **Setup** screen by selecting **User: Settings** from the main menu.
3. Click the **Add** button and enter your network code.

Note: The network code is in the documentation included with the purchase of your TPTF product. The network code is different from the main code. Once the code is accepted, your network will be enabled to run multiple concurrent copies.

* If you do not see the Network Client checkbox on your install make sure you have the latest install. The latest versions can be downloaded at www.TotalPTFitness.com. Download the demonstration version of your TPTF product to get the latest install.

For Citrix, SQL, or terminal server installations, contact us at support@totalptfitness.com